
# venome®

Filler

## Venome L Lidocaine

CE 0297

**C & V Technologies GmbH**  
Höyendörffstraße 20a 16761 Hennigsdorf,  
Niemchchina  
**S & V Technologies GmbH**  
Neuendorffstraße 20a 16761 Hennigsdorf,  
Germany  
Tel.: +49 3302-202-7900  
Fax: +49 3302-202-7905  
e-mail: [info@skin-vision.eu](mailto:info@skin-vision.eu)

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## Venome Filler Venome L Lidocaine

Medical device according to 93/42/EEC

**Description**  
Venome L Lidocaine is a sterile, apyrogenic clear gel made of highly purified crosslinked hyaluronic acid in physiological buffer solution. It is of biofermentative origin and free of animal traces. The gel contains 0.3% lidocaine.

**Composition**  
1 ml contains:  
• sodium hyaluronate 24 mg  
• lidocaine hydrochlorid 3 mg  
• phosphate buffer 1 ml q.s.

**Content**  
Per sterile disposable syringe of Venome L Lidocaine: two sterile single use 27G 1/2 needles to be used exclusively for the injection of this Venome L Lidocaine; one Instruction for Use; a set of secondary labels for the patient and the patient's record.

**Indication**  
Viscoaugmentation for correction of volume losses or defects related to disease or age in the face, e.g. for volume augmentation of cheeks and for correction of medium to deep wrinkles like marionette lines and nasolabial folds. Lidocaine has been added to the gel to reduce the local sensation of pain after the injection.

### Mode of action

The injection of Venome L Lidocaine into the extracellular matrix adds volume to the tissue. Wrinkles, creases and folds are corrected. The effect may last between 6-9 months; however, it is influenced by the severity of volume loss or defect of the treated area, the technique of injection, as well as by biological, lifestyle, and environmental factors. Crosslinking of the hyaluronic acid slows down degradation.

### Application and Techniques

Intradermal injection of Venome L Lidocaine should be performed by physicians, with appropriate training in dermal filler injection. The correct injection technique is essential for the success of the treatment. Injection into the wrong skin layer, an inappropriate technique or a high injection rate may lead to pronounced side effects.

Venome L Lidocaine should be used at room temperature. It is suitable for the replacement or augmentation of tissue volume in order to compensate for volume loss caused by skin aging. It is intended for mid to deep dermal layer injection.

The volume to be injected per treated wrinkle differs individually and depends on the grade of the wrinkle. It is recommended to slightly undercorrect the wrinkle. Also in the case of wrinkles that have already been corrected once, consider the risk of overcorrection. Selection of the injection technique depends on the type and degree of the folds. Markedly deep folds that do not disappear by stretching will be difficult to lift up.

In order to prevent the formation of nodules, an even distribution of the gel should be ensured. Slow injection will help prevent or lessen injection related side effects.

Following the injection, the gel should be distributed in the tissue by gentle massage. It is recommended to cool the treated areas of the face for about 10 minutes following injection. After the initial treatment, touch-up injections may be necessary in some cases to achieve the desired correction of the wrinkles.

### Contraindications

Venome L Lidocaine is contraindicated in case of:  
• allergy to hyaluronic acid, lidocaine hydrochloride or local anesthetics of the amide type  
• auto-immune disease or during immunotherapy  
• acute skin infection in the area to be treated (e.g. acne)  
• acute influenza and feverish cold  
• acute systemic infection (herpes, psoriasis)  
• coagulation disorder  
• a tendency to form scars  
• pregnancy or breast-feeding  
• people under the age of 18

Venome L Lidocaine must not be injected into:

- orbita region, glabella and forehead
- blood vessels
- joints
- breasts

### Precautions

As a rule, the physician should carefully investigate the patient's history with respect to diseases, other anti-aging treatments and individual factors that may increase the risk of adverse reactions. The patient should be informed prior to treatment about the concept of tissue augmentation with hyaluronic acid, common injection-related reactions and any possible side effects.

There are no clinical data on the tolerability and effectiveness of Venome L Lidocaine in patients with pre-existing conditions, such as active autoimmune diseases and severe multiple allergies or in patients with a history of anaphylactic shocks.

The physician must evaluate each individual case to decide whether or not it is safe to inject Venome L Lidocaine, depending on the type of disease and the associated treatment. For these patients, it is recommended to propose an initial double test. Refrain from injection if the disease is active or evolving. Furthermore, it is recommended to carefully monitor these patients after the treatment.

Venome L Lidocaine should only be injected into healthy, non-inflamed and disinfected areas of skin. Wrinkles, creases and folds shall not be overcorrected. Excessive hyaluronic acid gel implantation may result in formation of nodules, large swellings and dislocation.

When using a sharp needle for gel implantation the fine blood vessels may be injured, inducing a hematoma. The inappropriate use of a blunt cannula may increase the trauma in the skin tissue, resulting in stronger swelling.

Patients should not take aspirin, ginkgo preparations or other blood thinning substances one week prior to treatment.

After treatment, the patient should not be exposed to excessive sun or intense heat (solarium, sun bath, sauna), at least until the initial reactions (redness, swelling) have resolved.

### Warnings

Venome L Lidocaine is only suitable for intradermal implantation. It must not be mixed with other filler products (e.g. permanent tissue fillers). Venome L Lidocaine should not be combined with other methods of aesthetic medicine, such as peeling, dermabrasion or any form of treatment with radiation (laser, radiofrequency, ultrasound) while the effects lasts (6-9 months). Venome L Lidocaine should not be used together with products that contain quaternary ammonium compounds because of possible incompatibilities. Before use, check the integrity of the packaging and the expiry date of the product. Do not use when packaging is compromised. If the 27G 1/2 needle gets blocked, the injection should be interrupted and the needle should be replaced.

For single use only. Discard any unused amounts. Re-use may pose a risk of infection caused by non-sterile product. Do not resterilize. The packaging is not autoclavable. Do not use in the operating room, because the outside of the syringe is not sterile.

Do not add other substances to Venome L Lidocaine, and do not transfer the product from the syringe into another container.

### General information on possible lidocaine side effects and restrictions for use

The total dose of administered lidocaine should be considered for concurrent treatments of different formulations (topical application). This also applies to simultaneous use with other local anaesthetics or substances with structural similarities to local anaesthetics of the amide type. The application may also result in a cumulative systemic toxicity.

The maximum dose per application of lidocaine HCl (without epinephrine) should be considered in case of a pretreatment with lidocaine. 300 mg (or 4.5 mg / kg) should not be exceeded in normal healthy adults. Overdose of lidocaine HCl regularly leads to symptoms in the central nervous system or to cardiovascular toxicity. During concurrent acidosis, hyperkalaemia, hypercalcaemia and hypoxia the toxic effects of local anaesthetics may be intensified and prolonged. Upon intoxication with local anaesthetics, centrally acting anaesthetics are contraindicated.

The physician must carefully evaluate each case of epilepsy, heart disease, severe hepatic impairment, severe renal impairment or porphyria, in order to determine whether or not treatment with Venome L Lidocaine is safe. Hereby, the severity of the disease and the type of the associated treatment must be considered.

Special caution is required for patients with congenital methaemoglobinemia, glucose-6-phosphate dehydrogenase deficiencies and for patients taking methaemoglobin-inducing drugs.

For athletes, it has to be noted that lidocaine treatment may cause a positive result in an anti-doping test.

### Side effects

Immediately after treatment, injection related reactions like redness and swelling and in rare cases itching or bruising may be observed at the injection site(s). The day after, a pronounced swelling may occur but this will resolve spontaneously within 3-7 days. Allergic reactions in the form of urticaria, oedema, bronchospasm, anaphylactic reactions/shocks (i.e. acute respiratory distress syndrome, cardiovascular system reactions) are extremely rare, because Venome L Lidocaine is protein free due to the non-animal source.

Rarely, formation of nodules or oedema may occur at the injection site.

In case of an inflammatory reaction it is recommended that the patient visits the physician.

Literature on hyaluronic acid based filler materials reports individual cases of granuloma, abscesses, delayed hypersensitivity reactions and discolouration. Other rare side effects found in the literature are necrosis caused by vessel occlusion. These are caused by the injection of the filler material into vessels causing occlusion, or when the filler material presses against the arterial wall.

Injection of lidocaine may cause a local redness or hypersensitivity reaction.

### Storage

Between 2°C and 25°C, in a dry place protected from sun light.

### Disposal

Directly after use dispose of syringe with unused gel and needles.

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