

Composition

- **Hyaluronic acid** - high molecular weight,
- **Rh collagen** - induction of procollagen synthesis,
- **Carboxymethylcellulose** - anti-aging effect, inhibition of the production of natural hyaluronidase.

The Rh collagen contained in the product induces the synthesis of procollagen in fibroblasts, thanks to which the skin regains elasticity.

Hyaluronic acid with a high molecular weight intensively moisturizes the skin and facilitates the proliferation of fibroblasts.

Carboxymethylcellulose has anti-aging properties and inhibits hyaluronidase, preventing the loss of natural hyaluronic acid.

SPECIAL WARNINGS/PRECAUTIONS FOR USE

- Use the standard precautions required for intradermal infiltration.
- As with any transdermal procedure, use of the device carries a risk of infection.
- A slight burning sensation may occur during treatment, which subsides within a few minutes. Local reactions usually associated with injection may occur in the treated area, such as swelling and erythema, which resolve spontaneously within a few hours after the procedure.
- For the success of the procedure, the physician's technical skills are essential.
- Any modification or use of the device, other than as specified in these instructions, may result in loss of sterility and its effectiveness. In such situations, the sterility guarantee the sterility and effectiveness of the device.
- Before starting the procedure, inform the patient about the indications and contraindications of the device, its incompatibilities and any adverse effects that may occur.
- Before starting the injection, thoroughly disinfect the area to be treated treatment. Inject (administer) slowly. Failure to do so may result in the risk of needle leakage and/or loss of the product at the Luer-Lock type connector.
- The amount of product to be injected depends on the type of area to be corrected.
- After injection, massage the area in question to ensure even distribution of the product.

CONTRAINDICATIONS

- DO NOT USE THE PRODUCT FOR APPLICATIONS OTHER THAN THOSE LISTED IN THESE INSTRUCTIONS.
- Do not inject the KARISMA product into the periorbital (eyelids, crow's feet) and glabellar (forehead) areas.
- Do not inject into blood vessels (intravascular).
- Do not inject excessive amounts of the product at the correction site (do not overcorrect).
- Do not use in case of individual hypersensitivity to the device.
- Do not use in patients diagnosed with vascular collagenopathies and autoimmune connective tissue diseases (connective tissue diseases).
- Do not use the device if there is infection or inflammation of the skin at or near the injection site.
- Do not use the device during pregnancy and lactation, although no specific contraindications are known.
- Do not inject the device in patients suffering from epilepsy that is not controlled medication (not being treated).
- Do not inject the device in patients prone to hypertrophic scar formation
- Do not inject in patients with known hypersensitivity to substances contained in the device.
- Do not use in patients under 18 years of age.
- Do not use the product directly or in combination with laser treatment, deep chemical peeling or dermabrasion. For superficial peeling, injection is discouraged if severe inflammatory reaction has occurred.

WARNINGS

The contents of the ampoule-syringe are sterile. The syringe is packaged in a tightly sealed unit packaging of blister type. The outer surface of the syringe is not sterile. Do not use after the expiration date indicated on the package. Do not use if the package is open or damaged.

The injection site must be within healthy skin. Do not inject into blood vessels, muscles, tendons or for breast augmentation. Do not use the filler with other products/devices. Do not re-sterilize. The product is intended for single use only. Store at 8 to 25°C, away from heat sources. Do not freeze. After use, immediately dispose of in accordance with applicable regulations. Keep out of reach of children.