

Advantages:

- Less pain and discomfort for the patient
- Minimal bleeding
- Fewer injection marks
- Easy to use – the blue line pointing upwards indicates the direction of the injection
- Optimal shape – can be placed on any standard syringe
- Certified medical device

CONTRAINDICATIONS

- The NANOSOFT device should not be used on broken skin, open wounds, cuts, or scars.
- The NANOSOFT device should not be used on rashes, skin infections, or any area where the skin is broken or altered.

UNWANTED RESULTS

- The observed local side effects, including swelling, redness and discoloration of the skin at the injection site, are most likely related to the injected substance and are not caused by the application of the NANOSOFT device.
- Adverse events associated with the use of the NANOSOFT device include a small risk of micro-bleeding.
- Other possible adverse events associated with the device may include, in rare cases, local infection at the injection site.