

Contraindications

It should not be used in the following cases:

- Allergy to hyaluronic acid, lidocaine hydrochloride, or other amide-type local anesthetics;
- Autoimmune diseases or during immunotherapy;
- Acute skin inflammation at the planned injection site (e.g., acne);
- Influenza and acute colds with fever;
- Acute systemic infection (herpes, psoriasis);
- Blood clotting disorders;
- Tendency to scarring;
- Pregnancy and breastfeeding;
- Under 18 years of age;

Venome Filler Lido must not be injected into:

- Blood vessels,
- Joints,
- Breasts.

The product is not recommended for use:

For treatments in the orbital area and transverse forehead wrinkles.

